DH-Ge detector (incl. -cal, -QC and -LP options) User notes



Overview

The Bentham DH-Ge houses a 2mm diameter active region Germanium detector, which may either mounted to the exit slit of Bentham monochromators, or used free- standing.

This device is responsive over the spectral range 800-1800nm.

The DH-Ge can be operated in either the DC or AC regimes; in the latter instance, a chopping frequency of 225Hz is suggested.



Should a calibration be purchased (_cal), please see the accompanying certificate for calibration results.

Electrical

This device is operated in the photovoltaic mode, without bias, and may be input into any Bentham current amplifier/ ORM via BNC. The pin configuration is inner pin cathode, outer sheath anode.

Specifications

The device specifications are as follow:-

Active area	Rec. range µm	Peak Resp. (A/W)	Shunt resistance kΩ	D* at lp cm.Hz ^{1/2} /W
2mm dia.	0.8 - 1.8	0.90	>90	5x10 ¹²

WEEE statement:

Bentham are fully WEEE compliant, registration number is WEE/CB0003ZR. Should you need to dispose of our equipment please telephone 0113 385 4352 or 4356, quoting account number 135419.



Bentham Instruments Limited 2 Boulton Road, Reading, Berkshire, RG2 0NH, UK Tel: +44 (0)118 975 1355 Fax: +44 (0)118 931 2971 Email: sales@bentham.co.uk Internet: www.bentham.co.uk