TEL301 Telescope User notes



Overview

The TEL301 is a direct view telescope, used to provide a welldefined field of view for the measurement of radiance and radiant intensity.

This device images the measurement plane onto an aperture; light incident on the aperture is output for measurement, the remainder of the view is reflected up to an eye-piece. In this manner, the user views a scene with a black spot corresponding to the area being measured.

The solid angle of measurement is defined by the diameter of the entrance lens (50mm) and the distance to the source.

Mechanical

The TEL301 may be mounted directly to any Bentham monochromator via an adaptor plate, or light may be transported via a fibre bundle which adapts to the telescope output. A 1/4" threaded hole is provided to the base of the telescope for mounting on a tripod.

Lenses

The following lenses are available for use with the TEL301.

| Lens | Wavelength range | Working distance |
|----------|------------------|------------------|
| TL1 | 300nm - 3µm | 1m - infinity |
| TL1(ACH) | 350nm - 1200nm | 1m - infinity |
| TL1(Q) | 200nm - 3µm | 1m - infinity |
| TLCaF2 | 200nm - 10µm | 1m - infinity |
| TL2 | 300nm - 3µm | 150mm |
| TL2(ACH) | 350nm - 1200nm | 150mm |
| TL2(Q) | 200nm - 3µm | 150mm |
| TL3 | Visible | 15mm |
| TL-MACRO | 300 - 3µm | 25mm |

All but the macro lens is mounted in a lens barrel on the end of a helical to adjust the lens position to achieve focus on the measurement plane. The position of the lens in the lens barrel can be adjusted depending on the working range used; closer to the aperture working towards the longer distances, further away from the aperture for closer work.

Field of View

A selection of four interchangeable apertures are provided to define the measurement field of view. These provide FOVs of ~6', 20', 1° and 2°, however the actual field of view (full angle) in radians is equal to the ratio of the aperture diameter to the distance between entrance lens and aperture.

The TEL301 provides linear coverage at a distance D of 2d tan (FOV/2).

The apertures are simply installed, with the inclined face of the aperture facing the source under measurement.

WEEE statement:

Bentham are fully WEEE compliant, registration number is WEE/CB0003ZR. Should you need to dispose of our equipment please telephone 0113 385 4352 or 4356, quoting account number 135419.



Bentham Instruments Limited 2 Boulton Road, Reading, Berkshire, RG2 0NH, UK Tel: +44 (0)118 975 1355 Fax: +44 (0)118 931 2971 Email: sales@bentham.co.uk Internet: www.bentham.co.uk