

# SPECTRORADIOMETER SOFTWARE DEVELOPMENT KIT

Release 3.1.01.0

Bentham Instruments Ltd 2 Boulton Road, Reading, Berkshire, RG2 0NH Tel: +44 (0)118 975 1355, fax: +44 (0)118 931 2971, email: sales@bentham.co.uk

Bentham Instruments Spectroradiometer Software Development Kit Version 3.1.01.0

First printed January 1995, 1999, this revision Jun 2011

Copyright © 2006-2011 by Bentham Instruments Ltd. All rights reserved

Purchasers of this product may copy and use the programming examples included on the disk, and other files where indicated. Purchasers may also make one copy of the disk for backup purposes. The software may not be copied or distributed in any other way.

Purchasers of this product may copy and use the programming examples included in this manual. No other parts of this manual may be reproduced or transmitted in any form or by any means, electronic, optical or mechanical, including photocopying and recording, or by any information storage and retrieval system, without permission in writing from Bentham Instruments Ltd.

# CONTENTS

INTRODUCTION	1
1. GETTING STARTED	3
1.1.       INDEX OF FILES ON THE DISK         1.2.       INSTALLATION         1.3.       MINIMUM SYSTEM REQUIREMENTS	4
2. THE HARDWARE CONTROL DLL	5
<ul> <li>2.1. WHAT IS A DLL</li> <li>2.2. THE SYSTEM MODEL</li> <li>2.2.1. List of Component types</li> </ul>	5 8
<ul><li>2.2.2. Component Types Alphabetical List</li><li>2.3. COMPONENT GROUPS</li><li>2.4. DLL FUNCTIONS</li></ul>	<i>10</i> 15
<ul> <li>2.5. HARDWARE ATTRIBUTE TOKENS</li> <li>2.6. DLL ERROR CODES</li> <li>2.7. CONTROLLING HARDWARE VIA THE DLL</li></ul>	44
<ul> <li>2.7.1. M300/Mc300 and DM150/DMc150 Monochromators</li> <li>2.7.2. TM300/TMc300, DTM300/DTMc300 and TTM300 Monochromators</li> <li>2.7.3. Filter Wheel</li> </ul>	45 46
2.7.4. SAM 2.7.5. TLS 2.7.6. MVSS	47
<ul> <li>2.7.7. SOB</li> <li>2.7.8. 228, 228A, 485 and 487 ADCs</li> <li>2.7.9. 225, 265, 267, 277, 485,487 and 477 Amplifiers</li> </ul>	48
<ul> <li>2.7.10. System Attributes</li> <li>2.8. USING THE DLL</li></ul>	51
3.1.       WINDOWS ERRORS         3.2.       HARDWARE CONTROL PROBLEMS	53

# INTRODUCTION

The Bentham Instruments Spectroradiometer Software Development Kit (SDK) is aimed at those who wish to develop their own applications to work with Bentham Instruments light measurement systems. The SDK centres on the use of the hardware control dynamic link library (DLL). The DLL provides a small but powerful suite of high-level functions that allow the developer to forget the low-level complexities of operating their hardware and concentrate instead on data acquisition and manipulation. The DLL is built from the code used in Bentham's own spectroradiometer control packages, meaning that it has been extensively tried and tested in a huge variety of systems.

This manual describes the use of the SDK. It is organised as follows:

Chapter 1 outlines the contents of the SDK and gives guidelines for installation.

**Chapter 2** describes the spectrometer control DLL in detail, including services that it offers and how they are used.

**Chapter 3** describes problems that you may encounter using the Spectroradiometer Control DLL and how to solve them.

#### 1.1. INDEX OF FILES ON THE DISK

#### \code examples \C

\Delphi \Java \LabView **MATLAB** \Python \VBA

Source code examples for the above languauges.

#### \hardware

1.

system.cfg	-	system configuration file (where provided)
system.atr	-	system configuration file (where provided)

#### \lib\

benhw32_fastcall.dll -	32 bit h	nardware control DLL
benhw32_stdcall.dll	-	<ul> <li>( Delphi _fastcall calling convention)</li> <li>32 bit hardware control DLL</li> <li>( Windows _stdcall calling convention)</li> </ul>
benhw32_cdecl.dll	-	32 bit hardware control DLL
ieee_32m.dll	-	( C/C++ _cdecl calling convention) 32 bit PC488 driver (4 versions included)

#### \lang\c

bendll.h	-	C header file
dllerror.h	-	C error code file
dlltoken.h	-	C attribute identifier file

#### \lang\generic

00		
dllerror.txt	-	text error code file
dlltoken.txt	-	text attribute identifier file

#### \lang\pascal

bendll.pas	-	Pascal DLL import unit
dllerror.pas	-	Pascal error code file
dlltoken.pas	-	Pascal attribute identifier file

#### \lang\vb

bendll.bas	-	VB DLL function declaration module
dllerror.bas	-	VB error code file
dlltoken.bas	-	VB attribute identifier file

#### 1.2. INSTALLATION

**Win2000/WinXP/Win7** - copy the contents of \lib bit on the SDK disk to your system directory (typically c:\windows\system), or the directory where your application resides.

If your hardware is not PC488 controlled system then the ieee\_32m.dll is a 'dummy' file but still needs to be installed. If you have a PC488 card then the copies and ieee\_32m.dll will be the same as those on the disks accompanying your card; licensing information can be found in the PC488 manual.

Where possible a system configuration file, \hardware\system.cfg, written specifically for your hardware is included. This file is essential for using the spectroradiometer control DLL and may be copied as required. Unless new hardware is added, this file should never be edited.

To use the example programs you should have installed the appropriate libraries as described above. The example programs also expect the system configuration file system.cfg to be in the same directory as they are launched from.

#### **1.3. MINIMUM SYSTEM REQUIREMENTS**

Intel Pentium II / 300 MHz or better Microsoft Windows 2000 Professional (SP4) or Windows XP Professional (SP2) 128Mb RAM or better CD-ROM or DVD-ROM Drive SVGA or higher-resolution monitor (XGA recommended)

# 2. THE HARDWARE CONTROL DLL

#### 2.1. WHAT IS A DLL

A dynamic link library (DLL) is an executable module that contains code or resources that can be used by clients such as other applications or DLLs. When a client wishes to utilise code or resources in a DLL, the DLL is loaded into memory and linked to the client at runtime. The DLL does not have to be written in the same language as the client and can be easily replaced with updated versions when necessary. For this manual a DLL will be considered as a precompiled library of functions that can be called by other applications.

#### 2.2. THE SYSTEM MODEL

In order for the DLL to be able to control a system it has to know what hardware it consists of. The DLL achieves this by constructing a **system model** that mirrors the hardware components. The system model is entirely modular, in the same way as the spectroradiometer system. It is the system model that co-ordinates the activities of the various pieces of hardware and allows almost any system, however large and complicated, to be controlled by the same set of simple, high-level functions.

The system model is built from a **system configuration** file. This describes the spectroradiometer in terms of what components it consists of, how the PC can communicate with them and how they interact. System configuration files rarely need to be changed, and where possible they are pre-written and supplied with the SDK.

Figure 1 shows a system configuration file that describes a simple TM300 based DC system.

#		
#Demo TM30	) system	
#		
USB COMMS	comms	vid 1240, pid 5892
USB	msd	vid 1240, pid 5893
Amp487	dc_amp	address 64
ADC487	adc	address 72
FW252	fwheel	drive 2
SAM	exit_sam	drive 1
TM300	mono	use fwheel, use exit_sam

Figure 1: Example System Configuration File

Figure 1 describes the following system:

1 USB\_COMMS - USB electronics controller at vid 1240 and pid 5892

1 USB – USB mechanical controller at vid 1240 and pid 5892

1 Amp487 - DC amplifier at address 64 on the I2C bus in the electronics bin

1 ADC487 - ADC at address 29 on the I2C bus in the electronics bin

1 252 filter wheel controlled by drive 2 of the MSD

1 SAM controlled by drive 1 of the MSD

1 TM300 monochromator housing the previously described SAM and filter wheel and MSD board

The format of the file is as follows:

- One component is described per line,
- For each line:
- The first column declares the type of the component,
- The second column gives the component an identifier,

- The third column contains any **parameters** (e.g. address), with multiple parameters separated by commas,

• Anything following a '#' on a line is a comment and is ignored.

The system configuration file is the means by which the DLL knows what hardware is present and how it is connected to the PC. It must contain one line for each PC controlled hardware module.

The first task when writing a system configuration is to declare how the PC communicates with any electronics in the system. Line 4 in figure 1 specifies that a USB\_COMMS controller with vid 1240 and pid 5892 is to be used (other options include a PC488 IEEE controller, a TAS016 RS232<>IEEE converter, INES PCMCIA IEEE controller or an RS232, RS232 controller).

The next task is to declare how the PC communicates with any mechanical components in the system. Line 5 in figure 1 specifies that a USB controller with vid 1240 and pid 5893 is to be used (other options include a PMC, MSD or MAC which can be RS232 or IEEE).

These lines defining how the PC communicates must appear in the file before any other system component so the DLL knows how to communicate with the hardware in the system.

The next step is to add a line describing each hardware module in turn. A good 'bottom-up' approach is to start with each electronics module (e.g.amplifiers and ADC), then each mechanical/stepper motor drive controlled unit (e.g. filter wheel, SAMs, SOBs, MVSSs) and finally the monochromator.

One component is described per line. The first column is the **type**. This tells the DLL what the component is, e.g. 487 amplifier, 487 ADC, PMC. The next section gives a complete list of component types. The type is case sensitive.

The second column is an **identifier** by which the component can be referred to. The identifier can contain any alphanumeric characters and the '\_' (underscore) character, but the first character must be a letter or '\_'. Each identifier must be unique and cannot be any of the system configuration keywords (i.e. a type or parameter). The identifier is case sensitive.

The third column is a list of component **parameters**. These provide details such as the address of the component; the next section gives the valid parameters for each component type. If there is more than one parameter they must be separated by commas. The parameter list is case sensitive.

In systems with more than one electronics controller the entry for each device must specify which controller the component is connected to. This is achieved by the use parameter. Inserting the identifier of a controller as the value for the use parameter instructs the system model to operate the component via that controller. For example, the configuration file for a system with an USB\_COMMS controlled 477Amp could contain:

USB\_COMMS comms2 vid 1240, pid 5891 477Amp pre amp address 104, use comms2

In systems with more than one stepper motor drive (SMD) the entry for each SMD controlled component must specify which SMD the component is connected to. This in a similar way to the electronics. For example, the configuration file for a system with an MSC1 controlled MVSS could contain:

MSC1	slit drive	address 20
MVSS	entrance slit	motor 1, use slit drive

This also specifies that the MVSS is to be addressed as motor 1. When there is more than one SMD in a system and the SMD is not specified for a component the DLL will assign a default, but it may not be the correct one.

The *use* parameter also allows some components (specifically filter wheels, SAMs and MVSSs) to be attached to the monochromator. The system model needs to know what components are part of the monochromator so that it can correctly co-ordinate all of their operations. This is illustrated in the final line of figure 1 which describes a TM300 that contains a SAM and filter wheel.

The use parameter is the reason for the suggested bottom-up order of adding components to the system configuration. This is because one component cannot be used by another unless it has been declared first.

Figure 2 is a more sophisticated example of a system configuration describing an AC system based on a DTM300 with a filter wheel, 2 SAMs and 3 MVSSs. The number of component types along with all of their different parameters may seem daunting, but once set up the system configuration file should never need altering. As the DLL functions access components via their identifiers as defined in the configuration file the only usual reason for looking at it is to find out what these identifiers are.

Once a system configuration file has been written it is passed to the DLL using the *Bl\_build\_system\_model* function. This instructs the DLL to compile a system model from the system configuration file. The function reports its success or failure and gives the reason for any failure to compile.

#\_\_\_\_\_\_ # Example system configuration file 2 # Specifies CEC488 controlled system PC488 comms address 21 # Defines an MSD3 connected to COM1 rather than the IEEE bus mono drive port COM1 MSD3 # Defines a MAC at address 20 slit drive address 20, cards 3 MAC # Defines a 277 pre-amplifier at address 28 Amp277 pre amp address 28 # Defines a 225 lock-in amplifier at address 27 Amp225 ac amp address 27 # Defines a 228A ADC at address 28 ADC228A adc address 29 # Defines 3 MVSSs, all controlled by the same MAC MVSSentrance\_slitdrive 1, use slit\_driveMVSSmiddle\_slitdrive 2, use slit\_driveMVSSexit\_slitdrive 3, use slit\_drive # Defines a filter wheel controlled by card 3 of the MSD FW252 fwheel card 3, use mono drive # Defines a SAM controlled as SAM1 on card 1 of the MSD exitl sam card 11, use mono drive SAM # Defines a SAM controlled as SAM2 on card 1 of the MSD exit2 sam card 12, use mono drive SAM # Defines a DTM300 containing all of the slits, SAMs and # filter wheel and driven by the MSD. Note that long lines are
# OK, but cannot contain any line breaks. DTM300 use fwheel, use exit1 sam, use exit2 sam, mono use entrance slit, use middle slit, use exit slit, use mono drive

Figure 2: Example System Configuration File

### 2.2.1. List of Component types

This section gives a complete list of component types. Each type is described, along with its parameters.

Component Type	Component Modelled	Parameters
PC488,	PC488 IEEE controller	address 030
CEC488		port \$ <i>address</i>
TAS016, RS232	TAS016 RS232<>IEEE converter	port COM <i>n</i>
USB_COMMS	USB Interface	pid <i>n (5892)</i> vid <i>n (1240)</i>
MSD	MSD micro-stepping drive	address 030
MOD		port COM <i>n</i> cards <i>n</i>
MAC	MAC micro-stepping drive	address 030 port COMn
USB	USB MAC micro-stepping drive	cards <i>n</i> pid <i>n (5893)</i> vid <i>n (1240)</i>
MSC1	MSC1 micro-stepping drive	address 030
PMC	PMC micr-stepping drive	address 030
ADC228	228 ADC	address 030
		use <i>identifier</i>
ADC228A	228A ADC	address 030 use <i>identifier</i>
ADC487	487 ADC	address 72,74,76,78
		use identifier
ADC485	485 ADC	address 96,98,100,102
		use identifier
Amp225	225 lock-in amplifier	address 030
		use identifier
Amp265	265 DC amplifier	address 030 use <i>identifier</i>
Amp267	267 DC amplifier	address 030
/////		use <i>identifier</i>
Amp277	277 pre-amplifier	address 030
		use identifier
Amp487	487 DC amplifier	address 64,66,68,70
Amp485	485 lock-in amplifier	use <i>identifier</i> address 88,90,92,94
		use identifier
Amp477	477 pre-amplifier	address
		104,106,108,110 use identifier
SB262	262 switch box	address 030
SB462	462 switch box	address 112118
FW252	252 filter wheel	card integer
		positions 012
		use identifier

#### Table 1: Complete List of Component Types

Anonymous IEEE device	address 030
	use identifier
Anonymous USB device	address 0238
	use identifier
Motorised stage	motor 13
	use identifier
MK1 Motorised Slit	resolution n
	drive 1max Drives
	use identifier
MK2 Motorised Slit	drive 1max Drives
	use identifier
MK3 Motorised Slit	drive 1max Drives
	use identifier
EBox Monitor	address 152
	use identifier
Swing Away Mirror	card integer
	use identifier
Switch-over Box	card integer
	use identifier
Triple Light Source	card integer
	use identifier
M300E monochromator	use identifier
	selfpark 1 or 0
M300HR monochromator	as M300E
M300T monochromator	as M300E
DM150 monochromator	as M300E
TM300 monochromator	use identifier
	subtractive 1 or 0
DTM300 monochromator	use identifier
	subtractive 1 or 0
TTM300 monochromator	use identifier
HR600 monochromator	use identifier
	subtractive 1 or 0
	MK1 Motorised Slit MK2 Motorised Slit MK3 Motorised Slit EBox Monitor Swing Away Mirror Switch-over Box Triple Light Source M300E monochromator M300HR monochromator M300HR monochromator DM150 monochromator TM300 monochromator

### 2.2.2. Component Types Alphabetical List

The following is an alphabetical list of all of the component types that can be used in the system model, along with valid parameters and any other notes.

ADC228	add a 2	228 ADC to the s	system model
Parame	eters	address use	the address of the 228 specify which controller to use
ADC228A	add a 2	228A ADC to the	system model
Parame	eters	address use	the address of the 228A specify which controller to use
ADC485	add a 4	485 ADC to the s	system model
Parame	eters	address use	the address of the 485ADC specify which controller to use
ADC487	add a 4	487 ADC to the s	system model
Parame	eters	address use	the address of the 487ADC specify which controller to use
Amp225	add a 2	225 AC amplifier	to the system model
Parame	eters	address use	the address of the amplifier specify which controller to use
Amp265	add a 2	265 DC amplifier	to the system model
Parame	eters	address use	the address of the amplifier specify which controller to use
Amp267add a 267 DC amplifier to the system model			
Parame	eters	address use	the address of the amplifier specify which controller to use
Amp277	add a 2	277 pre-amplifier	to the system model
Parame	eters	address use	the address of the amplifier specify which controller to use
Amp477	add a 4	477 pre-amplifier	to the system model
Parame	eters	address use	the address of the amplifier specify which controller to use
Amp485	add a 4	485 AC amplifier	to the system mod <b>el</b>
Parame	eters	address use	the address of the amplifier specify which controller to use
Amp487	add a 4	487 DC amplifier	to the system model
Parame	eters	address use	the address of the amplifier specify which controller to use

CEC488	add a	CEC488 (or com	patible) interface card
Param	eters	address port	the IEEE address of the card the memory address of the card (only required if card not at default address, e.g. INES PCMCIA IEEE card)
DM150	add a l	DM150 monochro	omator to the system model
Param	eters	use	add a slave component to the DM150 or specify
Notes		lfpark parameter	which SMD to use a value of 1 indicates self-park capability are filter wheels (5), SAMs (10) and MVSSs (5). is only used if the DM150 is PMC or MSC1
DTM300	add a l	DTM300 monoch	romator to the system model
Param	eters	use	adds a slave component to the DTM300 or specify which SMD to use
Notes	Valid s	lave components	are filters wheels (5), SAMs (10) and MVSSs (5).
EBox_Monitor	add ar	ebox monitor to	the system model
Param	eters	address use	the address of the ebox monitor specify which controller to use
FW252	define	a 252 filter whee	for use in a monochromator
Param	eters	use	the MAC, MSD or PMC that the filter wheel is connected to
Notes			which MSD card the filter wheel is connected to which MAC drive the filter wheel is connected to the number of filter positions (default is 6) a value of 1 indicates self-park capability is only used if it is PMC or MSC1 controlled. r is optional; the default (6) is used if omitted.
HR600	add a	HR600 monochro	omator to the system model
Param	eters	use	add a slave components to the HR600 or specify
Notes	Valid s	lave components	which SMD to use are filters wheels (5), SAMs (10) and MVSSs (5).
IEEEDevice	add ar	anonymous IEE	E Device
Param	eters	address use	the address of the device specify which controller to use
M300E	add ar	M300E monoch	romator to the system model
Param	eters	use	add a slave component to the M300E or specify which SMD to use
Notes		lfpark parameter	a value of 1 indicates self-park capability are filter wheels (5), SAMs (10) and MVSSs(5). is only used if the M300E is PMC or MSC1
M300HR	add an	M300HR monoc	hromator to the system model
Param	eters	use selfpark	add a slave component to the M300HR or specify which SMD to use a value of 1 indicates self-park capability
Notes	Valid s		are filter wheels (5), SAMs (10) and MVSSs(5).

The selfpark parameter is only used if the M300HR is PMC or MSC1 controlled. M300T add an M300T monochromator to the system model **Parameters** add a slave component to the M300T or specify use which SMD to use a value of 1 indicates self-park capability selfpark **Notes** Valid slave components are filter wheels (5), SAMs (10) and MVSSs(5). The selfpark parameter is only used if the M300T is PMC or MSC1 controlled. add a MAC to the system model MAC **Parameters** address address of the MAC cards number of motors the MAC can control MOTOR add a motorised stage to the system model **Parameters** which SMD the motor is driven by use which motor/drive the SMD refers to this as drive the number of steps per unit (e.g. steps/mm, steps/°) resolution MSC1 add an MSC1 to the system model **Parameters** address address of the MSC1 MSD add an MSD to the system model **Parameters** address of the MSD (when IEEE controlled) address serial port the MSD is connected to (when under port direct PC RS232 control) cards number of controller cards in the MSD **Notes** The address and port parameters are mutually exclusive; the MSD is either IEEE or RS232 controlled. MVSS-MK1..Mk3 define a motorised slit for use in a monochromator **Parameters** SMD which the slit is connected to use motor \_ which motor/drive of the SMD the MVSS is connected to Notes A motor drive must be specified with the use parameter. For a PMC controlled MVSS the motor parameter can be 1 or 2. For an MSC1 controlled MVSS it can be 1, 2 or 3. For MAC and IMAC controllers the motor parameter can be any of the available drives up to the max number of drives. PC488 add a PC488 (or compatible) interface card Parameters address the IEEE address of the card the memory address of the card (only required if port card not at default address, e.g. INES PCMCIA IEEE card) **PMC** adds a PMC to the system model **Parameters** address address of the PMC **RS232** communicate with hardware via RS232 **Parameters** port serial port to be used (e.g. COM1)

SAM	add a S	SAM to the syste	em model or define it for use in a monochromator
Param	eters	use	MAC, IMAC, MSD or PMC that the SAM is
		card drive	connected to which card the SAM is connected to (MSD control) which drive the SAM is connected to (MAC control)
Notes	drive/ca For exa the carc of opera	ard parameter sp ample, drive 21 r d, the second th	hat can control two SAMs or SOBs per card, the becifies the drive/card number and SAM number. refers to SAM1 on drive 2 of a MAC; the first digit is e SAM. For IMACs and MSDs that are only capable SOB per card, the card parameter is a single digit hber.
SB262	add a 2	er to the syster?	n model
Param	eters	address- use	address of the 262 specify which controller to use
SB462	add a 4	162 to the syster	n model
Param	eters	address use	address of the 462 specify which controller to use
SOB	add a S	SOB to the syste	m model or define it for use in a monochromator
Param		use card drive	MAC, IMAC, MSD or PMC that controls the SOB which card controls the SOB (MSD) which drive controls the SOB (MAC)
		unvo	
(card) ا drive 2 the SO	paramete 1 refers t B. For M	er specifies the o to SOB1 on driv /ISDs that are or	at can control two SAMs or SOBs per card) the drive drive (card) number and SOB number. For example, e 2 of a MAC; the first digit is the drive, the second hly capable of operating one SAM/SOB per card, the specifying the card number.
(card) ا drive 2 the SO	paramete 1 refers t B. For M arameter	er specifies the of to SOB1 on drive ASDs that are of is a single digit	at can control two SAMs or SOBs per card) the drive drive (card) number and SOB number. For example, e 2 of a MAC; the first digit is the drive, the second hly capable of operating one SAM/SOB per card, the
(card)   drive 2 the SO card pa	paramete 1 refers t B. For N arameter <i>commu</i>	er specifies the of to SOB1 on drive ASDs that are of is a single digit	at can control two SAMs or SOBs per card) the drive drive (card) number and SOB number. For example, e 2 of a MAC; the first digit is the drive, the second nly capable of operating one SAM/SOB per card, the specifying the card number.
(card)   drive 2 the SO card pa TAS016	paramete 1 refers t B. For M arameter <i>commu</i> eters	er specifies the of to SOB1 on drive ASDs that are of is a single digit unicate with hard port	at can control two SAMs or SOBs per card) the drive drive (card) number and SOB number. For example, e 2 of a MAC; the first digit is the drive, the second hly capable of operating one SAM/SOB per card, the specifying the card number.
(card)   drive 2 the SO card pa TAS016 Param	paramete 1 refers t B. For M arameter commu eters add a 1	er specifies the of to SOB1 on drive ASDs that are of is a single digit unicate with hard port	at can control two SAMs or SOBs per card) the drive drive (card) number and SOB number. For example, e 2 of a MAC; the first digit is the drive, the second hly capable of operating one SAM/SOB per card, the specifying the card number. <i>Iware via a TAS016</i> serial port that the TAS016 is connected to (e.g. COM1)
(card)   drive 2 the SO card pa TAS016 Param TLS	paramete 1 refers t B. For M arameter commu eters add a 7 eters	er specifies the of to SOB1 on drive ASDs that are on is a single digit <i>inicate with hard</i> port <i>Triple Light source</i> use card drive	at can control two SAMs or SOBs per card) the drive drive (card) number and SOB number. For example, e 2 of a MAC; the first digit is the drive, the second hly capable of operating one SAM/SOB per card, the specifying the card number. <i>Iware via a TAS016</i> serial port that the TAS016 is connected to (e.g. COM1) <i>Ce to the system model</i> SMD that the SAM is connected to which card the SAM is connected to (MSD control)
(card)   drive 2 the SO card pa TAS016 Param TLS Param	paramete 1 refers t B. For M arameter commu eters add a 7 eters add a 7	er specifies the of to SOB1 on drive ASDs that are on is a single digit <i>inicate with hard</i> port <i>Triple Light source</i> use card drive	at can control two SAMs or SOBs per card) the drive drive (card) number and SOB number. For example, e 2 of a MAC; the first digit is the drive, the second hly capable of operating one SAM/SOB per card, the specifying the card number. <i>Iware via a TAS016</i> serial port that the TAS016 is connected to (e.g. COM1) <i>ce to the system model</i> SMD that the SAM is connected to which card the SAM is connected to (MSD control) which drive the SAM is connected to (MAC control)
(card)   drive 2 the SO card pa TAS016 Param TLS Param TM300 Param	paramete 1 refers t B. For M arameter commu eters add a 7 eters add a 7 eters	er specifies the of to SOB1 on driv. ASDs that are of is a single digit unicate with hard port Friple Light source use card drive TM300 monochrous use	at can control two SAMs or SOBs per card) the drive drive (card) number and SOB number. For example, e 2 of a MAC; the first digit is the drive, the second hly capable of operating one SAM/SOB per card, the specifying the card number. <i>Iware via a TAS016</i> serial port that the TAS016 is connected to (e.g. COM1) <i>ce to the system model</i> SMD that the SAM is connected to which card the SAM is connected to (MSD control) which drive the SAM is connected to (MAC control) omator to the system model add a slave components to the TM300 or specify
(card)   drive 2 the SO card pa TAS016 Param TLS Param TM300 Param	paramete 1 refers t B. For M arameter <u>commu</u> eters <u>add a 7</u> eters <u>add a 7</u> eters Valid sl	er specifies the of to SOB1 on drive //SDs that are of is a single digit //is a single digit //inicate with hard port //iple Light source //is //inicate with hard port //inicate with hard //inicate with har	at can control two SAMs or SOBs per card) the drive drive (card) number and SOB number. For example, e 2 of a MAC; the first digit is the drive, the second hly capable of operating one SAM/SOB per card, the specifying the card number. <i>Iware via a TAS016</i> serial port that the TAS016 is connected to (e.g. COM1) <i>ce to the system model</i> SMD that the SAM is connected to which card the SAM is connected to (MSD control) which drive the SAM is connected to (MAC control) <i>omator to the system model</i> add a slave components to the TM300 or specify which SMD to use
(card)   drive 2 the SO card pa TAS016 Param TLS Param TM300 Param Notes	paramete 1 refers t B. For M arameter <u>commu</u> eters <u>add a 7</u> eters <u>add a 7</u> eters Valid sl <u>add a 7</u>	er specifies the of to SOB1 on drive //SDs that are of is a single digit //is a single digit //inicate with hard port //iple Light source //is //inicate with hard port //inicate with hard //inicate with har	at can control two SAMs or SOBs per card) the drive drive (card) number and SOB number. For example, e 2 of a MAC; the first digit is the drive, the second hly capable of operating one SAM/SOB per card, the specifying the card number. <i>Iware via a TAS016</i> serial port that the TAS016 is connected to (e.g. COM1) <i>ce to the system model</i> SMD that the SAM is connected to which card the SAM is connected to (MSD control) which drive the SAM is connected to (MAC control) <i>omator to the system model</i> add a slave components to the TM300 or specify which SMD to use s are filters wheels (5), SAMs (10) and MVSSs (5).

USB add USB mechanical controller to the system model

		pid vid	the product id of the USB hid device (default 5893) the vendor id of the USB hid device (default 1240)
USBDevice	add ai	n anony	mous USB device to the system model

Parameters	address	the address of the device
	use	specify which controller to use

**USB\_COMMS** add USB electronics controller to the system model

Parameters	pid	the product id of the USB hid device (default 5892)
	vid	the vendor id of the USB hid device (default 1240)

#### 2.3. COMPONENT GROUPS

Simple spectroradiometer systems contain the same basic set of components: a monochromator, a stepper motor drive, an amplifier, possibly a pre-amplifier and an ADC. Measurements are made by repeatedly selecting a wavelength and taking a reading. In this situation all of the components are used together.

More sophisticated systems are possible. For example, consider a system designed to take both AC and DC measurements, possibly during the same scan. This time measurements would be made by selecting a wavelength and taking a reading using just one of the sets of detection electronics. In this situation the system model needs to know what set of components it should be using to perform the required operation.

To allow this the system model includes the concept of **component groups**. A component group is a set of components that are used together to perform operations. At any given time there is one **active group**; this is the component group that is used for all operations, such as selecting a wavelength or taking a measurement. A component group is a sub-set of the system model.

New component groups are added with the *BI\_build\_group* function. This constructs a new, empty component group and returns its sequence number. The first group is group 1, the next group created is group 2, the next group 3, and so on. There can be up to 10 component groups. When the system model is built, the DLL constructs a single empty default component group (group 1).

Existing component groups can be edited by using the *BI\_group\_add* and *BI\_group\_remove* functions to add and remove components to and from groups. For each function the component is referred to by its identifier as set in the system configuration file and the group by its sequence number.

The function *BI\_use\_group* sets the active group. All of the hardware operation functions (*BI\_initialise*, *BI\_close\_shutter*, *BI\_park*, *BI\_zero\_calibration*, *BI\_select\_wavelength*, *BI\_autorange* and *BI\_measurement*) work with the active group. Obviously for simple systems the active group will always be the default group (group 1).

# Example sy	stem configuration			
PC488 MAC	comms smd	address address		cards 2
# Note DC an Amp267	d AC detection ele dc_amp	ctronics address	26	
Amp277 Amp225	pre_amp ac_amp	address address	-	
ADC228A	adc	address	29	
# SOB to sel SOB FW252		to ADC drive 11 drive 2		
TM300	mono	use fwhe	eel	

Figure 3: Example System Configuration Containing AC and DC Detectors

Figure 3 shows a system configuration for a TM300 based system with both AC and DC detection electronics. The input to the ADC (from the 267 or 225) is selected by the SOB,

which is controlled by the MAC. In order to make AC and DC measurements two groups, one containing the monochromator, MSD, SOB, ADC and 277 and 225, and one containing the monochromator, MSD, SOB, ADC and 267, must be created. This can be achieved using the following sequence of DLL function calls:

Assuming that when the SOB is relaxed/off the ADC is reading from the 225, the following sequence of commands is a very naive example of a simple measurement scan:

```
for wl = start wl to stop wl by increment
   begin
      {Measurement with group 1 (AC) }
      BI use group(1)
      BI_set( "adc_input", SOBState, 0, 0 )
      BI select wavelength ( wl, settle time )
      wait( settle time )
      BI autorange ()
      BI measurement ( reading )
      {Measurement with group 2 (DC) }
      BI use_group(2)
      BI set ( "adc input", SOBState, 0, 1 )
      BI select wavelength ( wl, settle time )
      wait( settle time )
      BI autorange()
      BI measurement ( reading )
   end
```

*BI\_zero\_calibration* needs to access the monochromator in order to correctly measure the zero offset and dark current for each state that the system will be in over the specified wavelength range. *BI\_measurement* also needs to access the monochromator in order to retrieve the zero calibration values for the current wavelength. This means that if measurements performed using a group are to be zero calibrated the group must include the monochromator.

While the DLL will not create default component groups, *BI\_save\_setup* records which components are in which groups in the system attributes file and *BI\_load\_setup* deletes any existing component groups and builds new ones from the information in the system attribute file.

#### 2.4. DLL FUNCTIONS

This section gives a complete description of each function in the spectroradiometer control DLL. The functions are listed in alphabetical order. C, Pascal and Visual Basic prototypes are given for each function. For those using different languages, the parameter types are as follows:

С	Pascal	Visual Basic	Description
int	integer	ByVal As Integer	16 bit signed integer
unsigned int	word	ByVal As Integer	16 bit unsigned integer
double	double	ByVal As Double	64 bit real number
char far *	Pchar	ByVal As String	pointer to a null terminated
			string
int far *	var int	As Integer	pointer to a 16 bit signed
			integer
unsigned int far *	var word	As Integer	pointer to a 16 bit unsigned
			value
long far *	var longint	As Long	pointer to a 32 bit signed
			integer

Table 1: DLL Function Parameter Types

For Visual Basic the symbol has been used as a continuation character. This indicates that a line of code has been split in the text but should be entered on one line in Visual Basic.

# BI\_automeasure

Syntax	BI_automeasure(reading)
Description	This function auto-ranges the amplifier(s) in the active group and returns the reading at the current wavelength. It takes the ADC offset and dark current (previously obtained by BI_zero_calibration) into account
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	reading - set to the reading at the current wavelength (if call is successful)
С	int far pascal BI_automeasure( double far *reading );
Pascal	function BI_ automeasure ( var reading : double ) : integer;
Visual Basic	Function BI_ automeasure ( reading As Double ) As Integer

# BI\_autorange

Syntax	BI_autorange
Description	This function auto-ranges the amplifier(s) in the active group.
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	None
С	int far pascal BI_autorange( void );
Pascal	function BI_autorange : integer;
Visual Basic	Function BI_autorange() As Integer

# BI\_build\_group

Syntax	BI_build_group
Description	This function constructs a new, empty component group. The DLL allows up to 10 component groups.
Return value	If the function was successful it returns the group number, otherwise it returns BI_error.
Parameters	None
С	int far pascal BI_build_group( void );
Pascal	function BI_build_group: integer;
Visual Basic	Function BI_build_group() As Integer

### Bl\_build\_system\_model

Syntax	BI_build_system_model( filename, error_report )
Description	This function instructs the DLL to build a system model from a specified system configuration file. This must be done before the DLL can perform any hardware control functions. The SDK is supplied with a customised system configuration file (system.cfg).
	This function <i>must</i> succeed before any other DLL function is called. If there is any error the function reports it via <i>BI_report_error</i> .
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	filename - the path name of the system configuration file error_report - contains error report if function failed (256 characters maximum)
С	int far pascal BI_build_system_model(char far *filename, char far *error_report );
Pascal	function Bl_build_system_model(filename: Pchar; error_report : Pchar) : integer;
Visual Basic	Function BI_build_system_model( ByVal filename As String, &ByVal error_report As String ) As Integer
Notes	In Visual Basic error_report must be declared as a fixed-length string, ie
	Dim error report As String * 256

# BI\_camera\_measurement

Syntax	BI_camera_measurement(id, wls, readings)		
Description	This function instructs the DLL to take a camera measurement using the camera defined by the id string. Two arrays are passed in, the first of which will be filled with the camera wavelengths, and the second the intensity readings. The camera measurement itself is performed by the external dll defined in the system configuration file for the relevant camera.		
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure		
Parameters	id- component identifier as set in the system configuration filewls- pointer to an array of doublesreadings- pointer to an array of doubles		
С	int far pascal BI_camera_measurement (char far *id double far *wls double far *reading);		
Pascal	function BI_ camera_measurement (id: pchar; wls: pdouble; readings: pdouble ): integer;		
Visual Basic	Function BI_ camera_measurement (ByVal id As String, wls As Double, readings As Double) As Integer		

### BI\_close

Syntax Bl\_close

**Description** This function instructs the DLL to destroy the system model and prepare for unloading.

Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	None
С	int far pascal BI_close( void );
Pascal	function BI_close: integer;
Visual Basic	Function BI_close () As Integer

# BI\_close\_shutter

Syntax	BI_close_shutter
Description	This function instructs the DLL to send the filter wheel in the monochromator of the active group to its shutter position.
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	None
С	int far pascal BI_close_shutter( void );
Pascal	function BI_close_shutter : integer;
Visual Basic	Function BI_close_shutter() As Integer

### BI\_component\_select\_wl

Syntax	BI_component_select_wl(id, wavelength, delay)
Description Return value	This function sends the specified component to the specified wavelength and recommends a settle delay time before any readings are taken. It does not perform the delay itself. The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	id - component identifier as set in the system configuration file wavelength - wavelength to go to (nm) delay - recommended delay (ms) before taking readings
С	int far pascal BI_component_select_wl (char far *id,double wavelength, long far*delay );
Pascal	function BI_component_select_wl (id: Pchar; wavelength: double; var delay: longint ) : integer;
Visual Basic	Function BI_component_select_wl (ByVal id As String, ByVal wavelength As Double, delay As Long ) As Integer

BI_delete_group	
Syntax	BI_delete_group(n)
Description	This function deletes the specified group from the system model
Return value	The return value indicates either the result of the call: BI_error - failure Or the number of groups remaining after the deletion.
Parameters	n - group number to be deleted
С	int far pascal BI_delete_group(int token)
Pascal	function BI_delete_group ( n : integer; ) : integer;
Visual Basic	Function BI_delete_group (ByVal n As Integer ) As Integer

BI_get	
Syntax	BI_get( id, token, index, value )
Description	This function returns the value of the specified component attribute.
Return value	The return value indicates the result of the call:BI_OK- successBI_invalid_token- invalid tokenBI_invalid_component-the component does not exist in this systemBI_invalid_attribute- the attribute is not accessible in this system
Parameters	<ul> <li>id - component identifier as set in the system configuration file</li> <li>token - the token for attribute to be retrieved</li> <li>index - for accessing indexed attributes and setups</li> <li>value - set to the value of the attribute (if call successful)</li> </ul>
С	int far pascal BI_get( char far *id int token, int index double far *value );
Pascal	function BI_get( id : Pchar; token : integer; index : integer; var value : double): integer;
Visual Basic	Function BI_get( ByVal id As String, ByVal token As Integer, &ByVal index As Integer, value As Double ) As Integer

BI_get_c_group	
Syntax	BI_get_c_group (group)
Description	This function returns the current group number
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	group - set to the current group index (if call is successful)
С	int far pascal BI_get_c_group ( int far * group);
Pascal	function BI_get_c_group ( var group: integer ) : integer;
Visual Basic	Function BI_get_c_group (group As integer ) As Integer

# BI\_get\_component\_list

BI_get_component_list (list)
This function returns a list of all components in the system model.
NOTE : this does not return components contained within monochromators, to retrieve this information run BI_get_mono_items on all monochromators in the system.
The return value indicates the result of the call: BI_OK - success BI_error - failure
list - a comma-delimited list of the identifiers of the components in the system model
int far pascal BI_get_component_list (char far *list);
function BI_get_component_list (list:pchar) : integer;
Function BI_get_component_list (ByVal list As String) As Integer
In Visual Basic description must be declared as a fixed-length string

BI_get_group	
Syntax	BI_get_group( group, description )
Description	This function returns a list of the identifiers of all of the components in the specified group.
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	group - the number of the group description - a comma-delimited list of the identifiers of the components in the group
С	int far pascal BI_get_group( int group, char far *description );
Pascal	function BI_get_group( group : integer; description : Pchar ) : integer;
Visual Basic	Function BI_get_group( ByVal group As Integer, ByVal description As String %) As Integer
Notes	In Visual Basic description must be declared as a fixed-length string

BI_get_hardware_type	
Syntax	BI_get_hardware_type (id, hardware_type)
Description	This function returns the hardware type of the specified component
	NOTE : The hardware_type codes can be found in the dlltokens file.
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	id: component identifier as set in the system configuration file hardware_type: an integer token indicating the hardware type
С	int far pascal BI_get_hardware_type(char far*id, int far* hardware_type);
Pascal	function BI_get_hardware_type(id: pchar; hardware_type : pinteger ): integer;
Visual Basic	Function BI_get_hardware_type(ByVal id As String, &ByVal hardware_type As Integer) As Integer

# BI\_get\_max\_bw

Syntax	BI_get_max_bw (group,start_wl,stop_wl,bandwidth)
Description	This function returns the maximum bandwidth that the slits may be set to across the given wavelength range
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	group - group number currently in use start_wl - the start wavelength stop_wl - the stop wavelength bandwidth - the return value of maximum bandwidth across the wavelength range
С	int far pascal BI_get_max_bw(int group, double start_wl, double stop_wl, double far* bandwidth);
Pascal	BI_get_max_bw(group: integer; start_wl,stop_wl : double; var bandwidth : double):integer;
Visual Basic	Function BI_get_max_bw(ByVal group As Integer, ByVal start_wl As Double, ByVal stop_wl As Double, ByVal bandwidth As Double) As Integer

# BI\_get\_min\_step

Syntax	BI_get_min_step (group,start_wl,stop_wl,min_step)
Description	This function returns the minimum step that the monochromator is capable of, given the wavelength range and grating used
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	group - group number currently in use start_wl - the start wavelength stop_wl - the stop wavelength min_step - the return value of min_step across the wavelength range (in nm)
С	int far pascal BI_get_min_step(int group, double start_wl, double stop_wl, double far* min_step);
Pascal	BI_get_min_step(group: integer; start_wl,stop_wl : double; var min_step : double):integer;
Visual Basic	Function BI_get_min_step(ByVal group As Integer, ByVal start_wl As Double, ByVal stop_wl As Double, ByVal min_step As Double) As Integer

# BI\_get\_mono\_items

Syntax	BI_get_mono_items (id, list)
Description	This function returns a list of all components contained within the specified monochromator
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	id: component identifier as set in the system configuration file list - a comma-delimited list of the identifiers of the components in the monochromator
С	int far pascal BI_get_mono_items(char far*id char far*list);
Pascal	function BI_get_mono_items (id: pchar; list: pchar): integer;
Visual Basic	Function BI_get_mono_items(ByVal id As String,
Notes	SyVal list As String) As Integer In Visual Basic description must be declared as a fixed-length string

# Bl\_get\_n\_groups

Syntax	BI_get_n_groups(n)
Description	This function returns the number of groups in the system model
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	n - number of groups
С	int far pascal BI_get_n_groups(int far* n);
Pascal	function BI_get_n_groups(var n: integer) : integer;
Visual Basic	Function BI_get_n_groups(ByVal n As Integer) As Integer

### BI\_get\_no\_of\_dark\_currents

Syntax	BI_get_no_of_dark_currents (num)
Description	This function returns the number of dark currents taken during the last zero calibrate
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	num: the number of dark currents taken during the last zero calibrate for the current group
С	int far pascal BI_get_no_of_dark_currents (int far* num);
Pascal	function BI_get_no_of_dark_currents (var num: integer): integer;
Visual Basic	Function BI_get_no_of_dark_currents (ByVal num As Integer) As Integer

# BI\_get\_str

Syntax	BI_get_str( id, token, index, s )
Description	This function returns the value of the specified component attribute, where the value is a string
Return value	The return value indicates the result of the call:BI_OK- successBI_invalid_token- invalid tokenBI_invalid_component-the component does not exist in this systemBI_invalid_attribute- the attribute is not accessible in this system
Parameters	<ul> <li>id - component identifier as set in the system configuration file</li> <li>token - the token for attribute to be retrieved</li> <li>index - for accessing indexed attributes and setups</li> <li>s - set to the string of the attribute (if call successful)</li> </ul>
С	int far pascal BI_get( char far *id int token, int index char far *s );
Pascal	function BI_get( token : integer; index : integer; s : Pchar) : integer;
Visual Basic	Function BI_get( ByVal id As String, ByVal token As Integer, & ByVal index As Integer, n As String ) As Integer

BI_get_zero	_calibration_info
Syntax	BI_get_zero_calibration_info (wavelength,darkcurrent,adc_offset)
Description	This function returns lists of wavelengths, dark currents and adc_offsets of all points from the last zero calibrate on the current group.
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	wavelength: pointer to an array of doubles that the function fills with the wavelengths of zero calibrations darkcurrent: pointer to an array of doubles that the function fills with the dark currents from zero calibration points adc_offset: pointer to an array of doubles that the function fills with the adc offsets from zero calibration points
С	int far pascal BI_get_zero_calibration_info (double far*wavelength, double far* darkcurrent, double far* adc_offset);
Pascal	function BI_get_zero_calibration_info (wavelength: pdouble darkcurrent: pdouble adc_offset: pdouble): integer;
Visual Basic	Function BI_get_zero_calibration_info (ByVal wavelength As Double, 失ByVal darkcurrent As Double

⇔ByVal adc\_offset As Double) As Integer

BI_group_add	
Syntax	BI_group_add( id, group )
Description	This function adds a component to the specified group. If the group already contains a component of the same type it is replaced by the new one.
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	<ul> <li>id - component identifier as set in the system configuration file</li> <li>group - number of the group to add the component to</li> </ul>
С	int far pascal BI_group_add(char far*id, integer group );
Pascal	function BI_group_add(id: Pchar; group: integer ) : integer;
Visual Basic	Function BI_group_add( ByVal id As String, ByVal group As Integer ) As Integer

# BI\_group\_remove

Syntax	BI_group_remove( id, group )
Description	This function removes a component from the specified group.
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	<ul> <li>id - component identifier as set in the system configuration file</li> <li>group - number of the group to remove the component from</li> </ul>
С	int far pascal BI_group_remove(char far*id, int group );
Pascal	function BI_group_remove(id: Pchar; group: integer ) : integer;
Visual Basic	Function BI_group_remove( ByVal id As String, ByVal group As Integer ) &As Integer

# **BI\_initialise**

Syntax	BI_initialise
Description	This function initialises the active group as follows:
SAMs 265 267	<ul> <li>set-up to operate monochromator</li> <li>sent to initial position</li> <li>set to start range</li> <li>set to start gain and channel for the current set-up</li> <li>set to start gain and channel for the current set-up</li> <li>set to start gain, channel, phase variable and quadrant, frequency mode and time constant for the current set-up</li> </ul>
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	None
С	int far pascal BI_initialise( void );
Pascal	function BI_initialise(): integer;
Visual Basic	Function BI_initialise() As Integer
NOTE :	This function <i>must</i> be successfully called before the DLL performs any other hardware control functions.

# BI\_load\_setup

Syntax	BI_load_setup( filename )
Description	This function sets up the DLL system model from a specified system attribute file (created using BI_save_setup). The SDK is supplied with an example system attribute file (system.atr).
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	filename - the path name of the system attribute file
С	int far pascal BI_load_setup( char far *filename );
Pascal	function BI_load_setup( filename : Pchar ) : integer;
Visual Basic	Function BI_load_setup( ByVal filename As String ) As Integer

# BI\_measurement

Syntax	BI_measurement( reading )
Description	Using the active group, this function returns the reading at the current wavelength. It takes the ADC offset and dark current (previously obtained by BI_zero_calibration) into account.
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	reading - set to the reading at the current wavelength (if call is successful)
С	int far pascal BI_measurement( double far *reading );
Pascal	function BI_measurement( var reading : double ) : integer;
Visual Basic	Function BI_measurement( reading As Double ) As Integer

# BI\_multi\_automeasure

Syntax	BI_multi_automeasure(reading)
Description	This function auto-ranges the amplifier(s) in each group and returns an array of readings at the current wavelength (one for each group). It takes the ADC offset and dark current (previously obtained by BI_multi_zero_calibration) into account
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	reading - set to the reading of the first group at the current wavelength (if call is successful). The reading for subsequent groups is retrieved by incrementing the pointer.
С	int far pascal BI_automeasure( double far* reading );
Pascal	function BI_ automeasure ( reading : pdouble ) : integer;
Visual Basic	Function BI_ automeasure ( reading As Double ) As Integer

# BI\_multi\_autorange

Syntax	BI_multi_autorange
Description	This function auto-ranges the amplifier(s) in all groups in the system.
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	None
-	None int far pascal BI_multi_autorange ( void );
C	
C Pascal	int far pascal BI_multi_autorange ( void );

# BI\_multi\_get\_no\_of\_dark\_currents

Syntax	BI_multi_get_no_of_dark_currents (group, num)
Description	This function returns the number of dark currents taken during the last zero calibrate
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	group- number of the group to get the number of dark currents from num: the number of dark currents taken during the last zero calibrate for the specified group
С	int far pascal BI_multi_get_no_of_dark_currents (int group,int far* num);
Pascal	function BI_multi_get_no_of_dark_currents (group : integer; var num: integer): integer;
Visual Basic	Function BI_multi_get_no_of_dark_currents (ByVal group As Integer &ByVal num As Integer) As Integer
BI_multi_get_zero_calibration_info	

Syntax	BI_multi_get_zero_calibration_info (group, wavelength, darkcurrent, adc_offset)
Description	This function returns lists of wavelengths, dark currents and adc_offsets of all points from the last zero calibrate on the specified group.
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	group- number of the group to get the zero calibration info from wavelength: pointer to an array of doubles that the function fills with the wavelengths of zero calibrations darkcurrent: pointer to an array of doubles that the function fills with the dark currents from zero calibration points adc_offset: pointer to an array of doubles that the function fills with the adc offsets from zero calibration points
C	int far pascal BI_multi_get_zero_calibration_info (int: group, double far*wavelength, double far* darkcurrent, double far* adc_offset);
Pascal	function BI_multi_get_zero_calibration_info (group : integer; wavelength: pdouble darkcurrent: pdouble adc_offset: pdouble): integer;
Visual Basic	Function BI_multi_get_zero_calibration_info (ByVal group As Integer &ByVal wavelength As Double, &ByVal darkcurrent As Double

SyVal adc\_offset As Double) As Integer

# BI\_multi\_initialise

DI_IIIulu_III		
Syntax Description	BI_multi_initialise This function initialises all group as follows:	
SAMs 265 267 277	- set to start gain and channel for the current set-up	
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure	
Parameters	None	
С	int far pascal BI_multi_initialise( void );	
Pascal	function BI_multi_initialise(): integer;	
Visual Basic	Function BI_multi_initialise() As Integer	
NOTE :	This function must be successfully called before the DLL performs any other multi group hardware control functions.	
BI_multi_measurement		
Syntax	BI_multi_measurement ( reading )	
Description	This function returns the readings at the current wavelength for all groups. It takes the ADC offset and dark current (previously obtained by BI_multi_zero_calibration) into account.	
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure	
Parameters	reading – pointer to an array of doubles in which the function stores the readings at the current wavelength (if call is successful)	
C	int far pascal BI_multi_measurement (double far *reading );	
Pascal	function BI_multi_measurement (reading : pdouble ) : integer;	
Visual Basic	Function BI_multi_measurement (reading As Double ) As Integer	

# BI\_multi\_park

Syntax	BI_multi_park
Description	This function parks the monochromator of all the groups if possible (e.g. a TM300, DTM300 or M300/DM150 capable of self-parking) and leaves the DLL and monochromator in a well-defined state. Where a monochromator can be parked this function must be called before BI_zero_calibration and BI_select_wavelength are used. BI_multi_park only needs to be called once; after this the DLL records the state of the monochromator(s).
	Parking a monochromator parks all turrets, the filter wheel and all MVSSs.
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	None
С	int far pascal BI_multi_park( void );
Pascal	function BI_multi_park : integer;
Visual Basic	Function BI_multi_park( ) As Integer

# BI\_multi\_select\_wavelength

Syntax	BI_multi_select_wavelength ( wavelength, delay )
Description	This function performs the following operations with all groups in the system:
	<ul> <li>sends the monochromator to the specified wavelength,</li> <li>selects the filter for the specified wavelength,</li> <li>positions all SAMs according to the specified wavelength,</li> <li>sets all MVSSs according to the specified wavelength</li> <li>configures all amplifiers for the specified wavelength</li> </ul>
	and recommends a settle delay time before any readings are taken. It does not perform the delay itself.
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	wavelength - wavelength to go to (nm) delay - recommended delay (ms) before taking readings
С	int far pascal BI_multi_select_wavelength(  double  wavelength, long far *delay );
Pascal	function BI_multi_select_wavelength( wavelength : double; var delay : longint ) : integer;
Visual Basic	Function BI_multi_select_wavelength( ByVal wavelength As Double, together the select_wavelength by the select select select and the select sel

## BI\_multi\_zero\_calibration

Syntax	BI_multi_zero_calibration ( start_wavelength, stop_wavelength )
Description	This function performs ADC offset and system dark current measurements for all group. The DLL records the values and uses them when calculating readings in BI_multi_measurement. To ensure accurate readings BI_multi_zero_calibration should be called before BI_multi_measurement is used; once before a wavelength range scan is normally sufficient.
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	start_wavelength – wavelength to zero calibrate from. stop_wavelength – wavelength to zero calibrate to.
С	int far pascal BI_multi_zero_calibration (double start_wavelength, double stop_wavelength );
Pascal	function BI_multi_zero_calibration (start_wavelength: double; stop_wavelength: double ) : integer;
Visual Basic	Function BI_multi_zero_calibration (ByVal start_wavelength As Double, &ByVal stop_wavelength As Double ) As Integer
BI_park	

Syntax	BI_park
Description	This function parks the monochromator of the active group if possible (e.g. a TM300, DTM300 or M300/DM150 capable of self-parking) and leaves the DLL and monochromator in a well-defined state. Where a monochromator can be parked this function must be called before BI_zero_calibration and BI_select_wavelength are used. BI_park only needs to be called once; after this the DLL records the state of the monochromator.
	Parking the monochromator parks all turrets, the filter wheel and all MVSSs.
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	None
С	int far pascal BI_park( void );
Pascal	function BI_park : integer;
Visual Basic	Function BI_park() As Integer

BI_read	
Syntax	BI_read( message, buffer_size, chars_read, id )
Description	This function reads message from the device at the specified IEEE address. It is intended to be used for communicating with devices that are not controlled by any other DLL functions (e.g. non-standard or third party hardware).
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	<ul> <li>message - a pointer to a buffer for storing any message from the device</li> <li>buffer_size - the maximum number of characters to read</li> <li>chars_read - the number of characters actually read</li> <li>id - component identifier as set in the system configuration file</li> </ul>
С	int far pascal BI_read( char far *message, word buffer_size, word far *chars_read, char far* id );
Pascal	function BI_read( message : pchar; buffer_size : word; var chars_read : word; id : pchar) : integer;
Visual Basic	Function BI_read( ByVal message As String, ByVal buffer_size As Integer,

Schars\_read As Integer, ByVal id As String ) As Integer

BI_report_error	
Syntax	BI_report_error
Description	This function returns an error code corresponding to the last hardware error encountered. Calling this function resets the error code in the DLL.
Return value	Indicates the last hardware error encountered.
Parameters	None
С	int far pascal BI_report_error( void );
Pascal	function BI_report_error : integer;
Visual Basic	Function BI_report_error() As Integer

## BI\_save\_setup

Syntax	BI_save_setup( filename )
Description	This function creates a system attribute file for the current state of the DLL system model. This file can be reloaded with BI_load_setup.
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	filename - the path name of the set up file
С	int far pascal BI_save_setup( char far *filename );
Pascal	function BI_save_setup( filename : Pchar ) : integer;
Visual Basic	Function BI_save_setup( ByVal filename As String ) As Integer

BI_select_w	vavelength
Syntax	BI_select_wavelength( wavelength, delay )
Description	Using the active group, this function performs the following operations:
	<ul> <li>sends the monochromator to the specified wavelength,</li> <li>selects the filter for the specified wavelength,</li> <li>positions all SAMs according to the specified wavelength,</li> <li>sets all MVSSs according to the specified wavelength</li> <li>configures all amplifiers for the specified wavelength</li> </ul>
	and recommends a settle delay time before any readings are taken. It does not perform the delay itself.
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	wavelength - wavelength to go to (nm) delay - recommended delay (ms) before taking readings
С	int far pascal BI_select_wavelength( double wavelength, long far *delay );
Pascal	function BI_select_wavelength( wavelength : double; var delay : longint ) : integer;
Visual Basic	Function BI_select_wavelength( ByVal wavelength As Double, toget & delay As Long ) As Integer

## BI\_send

BI_SCIIG	
Syntax	BI_send( message, id )
Description	This function sends a character string to the device at the specified IEEE address. It is intended to be used for communicating with devices that are not controlled by any other DLL functions (e.g. non-standard or third party hardware). It should <b>not</b> be used to communicate with devices that the DLL controls.
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	<ul><li>message - a pointer to the string to be sent to the device</li><li>id - component identifier as set in the system configuration file</li></ul>
C	int far pascal BI_send( char far *message, char far* id );
Pascal	function BI_send(message: pchar; id: pchar): integer;
Visual Basic	Function BI_send( ByVal message As String, ByVal id As String ) As Integer

BI_set	
Syntax	BI_set( id, token, index, value )
Description	This function sets the value of the specified component attribute.
Return value	The return value indicates the result of the call:BI_OK- successBI_invalid_token- invalid tokenBI_invalid_component-the component does not exist in this systemBI_invalid_attribute- the attribute is not accessible in this system
Parameters	<ul> <li>id - component identifier as set in the system configuration file</li> <li>token - the token for attribute to be set</li> <li>index - for accessing indexed attributes and setups</li> <li>value - the new attribute value</li> </ul>
С	int far pascal BI_set(char far*id, int token, int index, double value );
Pascal	function BI_set(id: Pchar;token: integer; index: integer; value: double ): integer;
Visual Basic	Function BI_set( ByVal id As String, ByVal token As Integer, ∜ByVal index As Integer, ByVal value As Double ) As Integer
Notes	Using <i>BI_set</i> never triggers a hardware settle delay. If an attribute is being changed that causes some hardware operation to occur is the responsibility of the client to ensure that any required settle delay is used.
BI_set_str	
Syntax	BI_set_str( id, token, index, s )
Description	This function sets the value of the specified component attribute, where value is a string
Return value	The return value indicates the result of the call:BI_OK- successBI_invalid_token- invalid tokenBI_invalid_component-the component does not exist in this systemBI_invalid_attribute- the attribute is not accessible in this system
Parameters	<ul> <li>id - component identifier as set in the system configuration file</li> <li>token - the token for attribute to be set</li> <li>index - for accessing indexed attributes and setups</li> <li>s - the new attribute string</li> </ul>
C	int far pascal BI_set(char far*id, int token, int index, char far* s );
Pascal	function BI_set(id: Pchar;token: integer; index: integer; s: Pchar ): integer;
Visual Basic	Function BI_set( ByVal id As String, ByVal token As Integer, ∜ByVal index As Integer, ByVal s As String ) As Integer

Bl_trace	
Syntax	BI_trace ( trace )
Description	This function toggles the dll tracer. When trace equals 1 the tracer will log dll function calls in C:\trace.txt
Return value	None
Parameters	trace: 0 – tracer off 1 – tracer on
С	void BI_trace ( int trace );
Pascal	procedure BI_trace(trace: integer);
Visual Basic	Sub BI_use_group( ByVal trace As Integer )

# BI\_use\_group

Syntax	BI_use_group( group )
Description	This function sets the active group.
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	group - the number of the group to make active
С	int far pascal BI_use_group( int group );
Pascal	function BI_use_group( group : integer ) : integer;
Visual Basic	Function BI_use_group( ByVal group As Integer ) As Integer

## **BI\_version**

Syntax	BI_version( s )
Description	This function returns version information for the DLL.
Return value	none
Parameters	s is a pointer to a buffer that the version information is to be copied to. The version string will be no longer than 80 characters.
С	void BI_version( char far *s );
Pascal	procedure BI_version( s : Pchar );
Visual Basic	Sub BI_version( ByVal s As String )

## BI\_zero\_calibration

Syntax	BI_zero_calibration( start_wavelength, stop_wavelength )
Description	This function performs ADC offset and system dark current measurements for the active group. The DLL records the values and uses them when calculating readings in BI_measurement. To ensure accurate readings BI_zero_calibration should be called before BI_measurement is used; once before a wavelength range scan is normally sufficient.
	The information is stored with the active group so it is possible to zero- calibrate each group and then swap between them at will; the DLL will always use the correct zero-calibration data.
Return value	The return value indicates the result of the call: BI_OK - success BI_error - failure
Parameters	start_wavelength – wavelength to zero calibrate from. stop_wavelength – wavelength to zero calibrate to.
С	int far pascal BI_zero_calibration( double start_wavelength, double stop_wavelength );
Pascal	function BI_zero_calibration( start_wavelength : double; stop_wavelength : double ) : integer;
Visual Basic	Function BI_zero_calibration( ByVal start_wavelength As Double, &ByVal stop_wavelength As Double ) As Integer

### 2.5. HARDWARE ATTRIBUTE TOKENS

Hardware attributes can be classified as one of two types:

i) **Structural attributes -** these describe what components the system is built from and how they interact. These attributes can be considered as hardwired.

ii) **Component attributes -** these describe how the components are set-up and used, which may differ from run to run.

The system configuration file passed to BI\_build\_system\_model contains structural attributes. It describes the hardware system in terms of its components and how they relate to each other. Examples of this type of information are:

- there is a 228A ADC at address 29 on the IEEE bus,
- there is a TM300 monochromator that contains a filter wheel controlled by card 2 in the MSD.

The file created by BI\_save\_setup consists of component attributes. It describes in detail exactly how each component is configured. Examples of this type of information are:

- the ADC is sampled 10 times per measurement,
- the recommended delay time after changing filter wheel position is 2 seconds.

Once *BI\_build\_system\_model* has been called the structural attributes cannot be changed. However, users may wish to alter component attributes such as ADC samples per reading or delay times. Component attributes are retrieved and set using the *BI\_get* and *BI\_set* functions.

To access an attribute with *BI\_get* or *BI\_set* the component id, attribute, and in some cases attribute index, must be specified. The component id is the same as that set in the system configuration file. The attribute token is a constant that the DLL recognises as referring to a particular attribute. The SDK includes a set of attribute token definition files for various languages. These give identifiers to the attribute tokens which should be used instead of the actual token values. This makes calls to *BI\_get* and *BI\_set* more readable and will prevent problems with future versions of the DLL where the token values may change. Table 2 gives a list of which languages are currently supported and how the attribute token definition files should be used.

The attribute index is used where an attribute token may refer to one of several values. For example, consider a 6 position filter wheel. The following function call will retrieve the filter value at position 3:

BI_get( "filter_wheel",	the identifier for the filter wheel
FWheelFilter,	attribute token for a filter value
З,	attribute index for position 3
filter_value )	attribute value returned in this variable

Language	Attribute Token Def File	Use
C++	dlltoken.h	#include "dlltoken.h" in any source file that uses BI_get or BI_set
Pascal	dlltoken.pas	Include dlltoken.pas in any source file that uses BI_get or BI_set
Visual	dlltoken.bas	Add dlltoken.bas in any project
Basic		that uses BI_get or BI_set

Table 2: Attribute Token Definition Files

If there is not a token definition file for a particular programming language the file dlltoken.txt contains a no-frills list of identifiers and values that can be converted to any language. This file may be copied as required.

The following is a list of all of the attribute tokens in the DLL. There is a brief description of each attribute including how it is indexed (if it is a multi-valued attribute), what values it can take and what it defaults to. The attributes are grouped by component.

Monochromator Attributes		
MonochromatorCurrentDialReading MonochromatorCurrentWL MonochromatorParkDialReading MonochromatorCurrentGrating MonochromatorScanDirection MonochromatorPark MonochromatorSelfPark MonochromatorModeSwitchNum MonochromatorModeSwitchState MonochromatorCanModeSwitch GratingA Gratingd GratingWLMax GratingWLMin GratingX1 GratingX2 GratingZ ChangerZ	current dial reading (M300/DM150) current wavelength (nm) park position (M300/DM150) current grating (1-3) scan direction repark monochromator legacy attribute (DM150) double – single sam index sam state boolean alpha (TM/DTM30) ruling density d (TM/DTM300) grating maximum $\lambda$ (TM/DTM300) grating constant 1 (HR600) grating constant 2 (HR600) grating constant 3 (HR600) zord (TM/DTM300) changer zord	0 undefined 0 undefined 1 -1 0 1 1 0 from d from d from d undefined undefined 0
Filter Wheel Attributes		
FWheelCurrentPosition FWheelFilter FWheelPositions	current filter position filter value number of filter positions	undefined 0 undefined
SOB Attributes		
SOBInitialState SOBState	state after BI_initialise current state	0 undefined
TLS Attributes		
TLSCurrentPosition TLSPOS TLSPositionsCommand TLSSelectWavelength TLSWL	current light source light source to switch to switch light source select wavelength for the TLS wavelength to switch light source at	1 0 undefined undefined 0
262 Attributes		
biRelay biCurrentRelay	relay status current relay status	0 undefined
SAM Attributes		
SAMInitialState SAMCurrentState SAMState SAMSwitchWL SAMDeflectName SAMNoDeflectName	state after BI_initialise SAMs current state state at SAMSwitchWI(n) state change wavelength (nm) name for deflected SAM state Name for undeflected SAM state	0 undefined 0 undefined Deflect No Deflect

## Stepper SAM Attributes

Stepper SAM Attributes		
SSEnergisedSteps SSRelaxedSteps SSMaxSteps SSSpeed SSMoveCurrent SSIdleCurrent	steps to energised position steps to relaxed position maximum steps motor speed motor current on move motor current when idle	undefined undefined undefined undefined undefined undefined
MVSS Attributes		
MVSSConstantBandwidth MVSSConstantwidth MVSSCurrentWidth MVSSPosition	current constant Bandwidth (nm) current constant width (mm) current width of the mvss (mm) set slit position in monochromator ('entrance', 'exit' or 'middle')	undefined 010 undefined 'entrance'
MVSSSetWidth MVSSSlitMode MVSSSwitchWL MVSSWidth	move the slit to the specified width (mm) the current slit drive mode state change wavelength (nm) width at specified state	010 0 / 1 undefined 010
ADC Attributes		
ADCAdaptiveIntegration ADCSamplesPerReading	the adaptive integration mode number of samples integrated for one reading	0/1 10
ADCSamplePeriod ADCVolts ADCAuxVolts ADCAuxOffset ADCAuxInput ADCTimeConstant ADCXYThetaReading	sample period in milliseconds take a reading read from a 487 Aux Input get/set the Aux Offset select the 487 Aux Input for use with Stanford lock-in for use with Stanford lock-in	100 undefined undefined undefined undefined
General Amplifier Attribute	9S	
AmpChannel AmpCurrentChannel AmpCurrentRange AmpGain	input channel (225, 267, 277 only) current channel (225, 267, 277 only) current range gain value	1 undefined undefined
AmpMinRange AmpMaxRange	minimum range maximum range	1 7 (225) 6 (others)
AmpOverload AmpOverrideWI	overload flag set group readings to 1 from this wavelength (nm)	0
AmpStartRange AmpUseSetup	starting range (265, 267, 277 only) $\lambda$ (nm) to switch to set-up (225, 267, 277 only)	1 0

 AmpCurrentSetup
 only)

 get/set current setup being used
 0

## 225 Attributes

A225fMode	frequency mode	1
A225PhaseQuadrant	phase quadrant	1
A225PhaseVariable	phase variable	0
A225TargetRange	target range	1
A225TimeConstant	time constant	1

Motorised Stage Attributes			
MotorPosition MotorStop	get/set motor position tell motor to stop moving	0 undefined	
EBox Monitor Attributes			
EBoxReadHV EBoxReadTemp EBoxReadHVRaw EBoxReadTempRaw EBoxWait EBoxRepeats EBoxCountsAtTargetTemp EBoxGradientTemp EBoxTargetTemp EBoxCountsAtTargetHV EBoxGradientHV EBoxTargetHV	Read High Voltage Read Temperature Read High Voltage without converting Read Temperature without converting Delay between HV or Temp readings Repeats of HV and Temp ADC reading at target temp Gradient of ADC counts vs Temp Target Temperature ADC reading at target HV Gradient of ADC counts vs HV Target High Voltage	Undefined Undefined Undefined 5000 1 11990 1288.19 23 10178 -216.76 750	
Camera Attributes			
CameraIntegrationTime CameraBeta CameraPhi	Get/Set Integration Time Get/Set Camera Beta Get/Set Camera Phi	Undefined 0 0	
Miscellaneous Attributes			
biSettleDelay	recommended settle delay (ms) after		
biMin biMax biParkPos biInput	operation minimum allowable position (Motor) maximum allowable position (Motor) position after parking (M300, DM150, Motor only) input channel indexed by setup (262,	-1 -1 -1	
biCurrentInput	265, 267, 277) current input channel (262, 265, 267,		
biMoveWithWavelength	277) change component state with		
biHasSetupWindow	wavelength ask whether setup window exists for component		
biHasAdvancedWindow biDescriptor	returns product description string (eg		
biProductName	DTM300) returns product name (eg DTM300 Monochromator)		
biParkOffset	park offset		
System attributes			
Sys225_277Input SysStopCount SysDarkIIntegrationTime	which 225 input the 277 output goes to stop-count value for AC zero-calibration Integration time to be used when calculating dark levels	1 1.0	
Bentham Hardware Types			
BenInterface BenSAM	hardware_type token for interfaces hardware_type token for SAMs		

BenSlit	hardware_type token for MVSSs
BenFilterWheel	hardware_type token for filter wheels
BenADC	hardware_type token for ADCs
BenPREAMP	hardware_type token for preamplifiers
BenACAMP	hardware_type token for ac amplifiers
BenDCAMP	hardware_type token for dc amplifiers
BenPostAMP	hardware_type token for post amplifiers
BenRelayUnit	hardware_type token for relay units
BenMono	hardware_type token for monochromators
BenAnonDevice	hardware_type token for anonymous devices
BenCamera	hardware_type token for camera devices
BenDiodeArray	hardware_type token for diode array devices
BenORM	hardware_type token for ORM400 devices
BenEBoxMonitor	hardware_type token for EBoxMonitor device
BenUnknown	hardware_type token for unknown devices

The use of attribute tokens in controlling components is explained in detail further on.

### 2.6. DLL ERROR CODES

Functions in the DLL return error codes to indicate their success or failure and the cause of any failure. These are defined as follows:

BI_OK -		function call succeeded	
BI_error -		function call failed	
Bl_invalid_token -	-	the function was passed an invalid attribute token	
BI_invalid_component the function was passed a component identifier that			
		does not exist	
Bl_invalid_attribute -	-	the function was passed an attribute token referring to an	
		attribute that does not exist or is inaccessible	

Hardware operation functions return either BI\_OK or BI\_error. When a hardware operation function returns BI\_error, BI\_report\_error can be used to return a hardware error code that describes exactly what caused the error. The hardware error codes are as follows:

accombes exactly what baased t	
BI_no_error -	no error to report
BI_PMC_timeout -	PMC not responding
BI_MSD_timeout -	MSD not responding
BI_MSC_timeout -	MSC1 not responding
BI_MAC_timeout -	MAC not responding
BI_MAC_invalid_cmd -	error in communication with MAC
BI_225_dead	225 not responding
BI_265_dead -	265 not responding
BI_267_dead -	267 not responding
BI_277_dead -	277 not responding
BI_262_dead -	262 not responding
BI_ADC_read_error -	ADC not responding
BI_ADC_invalid_reading	could not obtain valid ADC reading
BI_AMP_invalid_channel	invalid amplifier channel
BI_AMP_invalid_wavelength	invalid amplifier wavelength
BI_SAM_invalid_wavelength	invalid SAM wavelength
BI_MVSS_invalid_width	invalid MVSS wavelength
Bl_turret_invalid_wavelength	-attempt to send monochromator beyond wavelength range
BI_turret_incorrect_pos	-error in communication with MAC

Calling *BI\_report\_error* clears the error code, i.e. subsequent calls will return *BI\_no\_error* until another hardware error occurs.

The SDK includes error code definition files. These should be used in applications that need to know what the return values of DLL functions signify. Languages currently supported are:

Language	Error Code Definition File	Use
С	dllerror.h	Use #include "dllerror.h" in any source file that uses error codes
Pascal	dllerror.pas	Include dllerror.pas in any source file that needs to know the error codes
Visual Basic	dllerror.bas	Add in any project that uses error codes
Dasic	Table & Freese Oads D	

Table 3: Error Code Definition Files

If there is not an error code definition file for a particular programming language the file dllerror.txt contains a no-frills list of identifiers and values that can be converted to any language. This file may be copied as required

### 2.7. CONTROLLING HARDWARE VIA THE DLL

This section describes how attribute tokens are used to control the behaviour of the components in a system.

#### 2.7.1. M300/Mc300 and DM150/DMc150 Monochromators

The DLL is designed to allow monochromator control at the highest level; the function *Bl\_select\_wavelength* tells the DLL which wavelength the monochromator should be at and the DLL does the rest, co-ordinating the operation of the gratings and any filter wheels or SAMs. Before this can happen the DLL needs to know exactly how the monochromator is set up.

- **MonochromatorCurrentDialReading:** the current dial reading. For non-parking M300 and DM150 monochromators this **must** be set before any calls to *BI\_select\_wavelength*. Valid values are 0..999.99.
- **MonochromatorCurrentWavelength:** the wavelength (nm) that the monochromator is at; read-only.
- **MonochromatorScanDirection:** controls when anti-backlash precautions are taken. A value of 1 indicates that all wavelengths should be approached from a shorter wavelength, while 0 indicates that all wavelengths should be approached from a longer wavelength. If the DLL ever needs to approach a wavelength in the 'wrong' direction it will first overshoot the required position and then go to it from the other side, preventing any backlash error. Therefore if an application intends to scan from short to long wavelengths (e.g. UV to IR), MonochromatorScanDirection should be set to 1. If an application intends to scan from long to short wavelengths (e.g. IR to UV), MonochromatorScanDirection should be set to 0. The default value is 1 (increasing wavelength).

**Gratingd:** the line density d (lines mm<sup>-1</sup>) of the grating.

- **GratingWLMax:** the maximum wavelength (nm) that the grating should be used for. Setting this to 0 causes the DLL to use a default value chosen by the ruling density.
- **GratingWLMin:** the minimum wavelength (nm) that the grating should be used for. Setting this to 0 causes the DLL to use a default value chosen by the ruling density.
- **biSettleDelay:** the recommended settle delay after selecting a new wavelength. *Bl\_select\_wavelength* will return this value if the monochromator has changed position and it is the longest delay time triggered.

MonochromatorPark: called with any value, the monochromator will be reparked.

**MonochromatorSelfPark:** legacy token for telling the benhw whether a Dial Reading is entered by the user when the monochromator park is called, or whether the monochromator parks itself.

### 2.7.2. TM300/TMc300, DTM300/DTMc300 and TTM300 Monochromators

As with the M300 and DM150 monochromators, the DLL needs to know the set-up for TM300 and DTM300 monochromators. For a TM300 gratings are indexed as 1, 2 or 3. For a DTM300 gratings on turret 1 are indexed as 11, 12 or 13, and on turret 2 as 21, 22 or 23.

- **MonochromatorCurrentWavelength:** the wavelength (nm) that the monochromator is at. This attribute is read-only. Attempting to set it *will not* operate the monochromator; use the DLL function *BI\_select\_wavelength*.
- **MonochromatorScanDirection:** controls when anti-backlash precautions are taken. A value of 1 indicates that all wavelengths should be approached from a shorter wavelength, while

0 indicates that all wavelengths should be approached from a longer wavelength. If the DLL ever needs to approach a wavelength in the 'wrong' direction it will first overshoot the required position and then go to it from the other side, preventing any backlash error. Therefore if an application intends to scan from short to long wavelengths (e.g. UV to IR), MonochromatorScanDirection should be set to 1. If an application intends to scan from long to short wavelengths (e.g. IR to UV), MonochromatorScanDirection should be set to 0. The default value is 1 (increasing wavelength).

**GratingA:** the alpha value for a grating.

**Gratingd:** the line density d (lines mm<sup>-1</sup>) for a grating.

GratingWLMax: the maximum wavelength (nm) that the grating should be used for.

GratingWLMin: the minimum wavelength (nm) that the grating should be used for.

GratingZ: the zero order (zord) value for a grating.

**biSettleDelay:** the recommended settle delay after selecting a new wavelength. *BI\_select\_wavelength* will return this value if the monochromator has changed position and it is the longest delay time triggered.

### 2.7.3. Filter Wheel

Filter wheels may have 6, 8, 10 or 12 positions; 6 is standard. The number of positions is fixed in the system configuration file. Each position can be assigned a filter value. If the filter wheel is attached to a monochromator then a call to *BI\_select\_wavelength* will send the filter wheel to the position with the highest filter value that is less than or equal to the specified wavelength.

The last position is always used as the shutter position, but this may also be assigned a filter value. *Bl\_close\_shutter* will send the filter wheel to this position. Positions with no filters should have a filter value of 0 (default). The filter values do not have to be assigned in any particular order, but normal practice is to start with the lowest filter and work up.

- **FWheelCurrentPosition:** the current filter position. This *must* be set for non-parking M300 and DM150 systems before calling *BI\_zero\_calibration*, *BI\_select\_wavelength* or *BI\_close\_shutter*. Read only for MAC and MSD controlled monochromators
- **FWheelFilter:** filter value, indexed by position (1..number of positions). This value is used to select a filter when *BI\_select\_wavelength* is called. When there is no filter at position *n* FWheelFilter(*n*) should be 0 (default).
- **FWheelPositions:** specifies the number of filter positions. The last (highest numbered) position is always used as the shutter. Read only.
- **biSettleDelay:** the recommended settle delay after a change of filter position. *BI\_select\_wavelength* will return this value if the filter wheel has been moved and it is the longest delay time triggered.

### 2.7.4. SAM

SAMs can be in one of two defined states, energised or relaxed. SAMs have an initial state that they go to when *Bl\_initialise* is called and also maintain a list of up to 10 records that specify what state the SAM should be in at a given wavelength. If a SAM is attached to a monochromator then a call to *Bl\_select\_wavelength* will cause the SAM to go to the state corresponding to the record with the highest wavelength that is less than or equal to the wavelength specified.

- **SAMInitialState:** the state of the SAM after calling *BI\_initialise*. 0 indicates relaxed and 1 is energised. The default is 0 (relaxed).
- **SAMSwitchWI, SAMState:** used to access the SAM's wavelength-state records. Both attributes are indexed by the setup number (1..10). *SAMSwitchWI* refers to the wavelength and *SAMState* to the corresponding state.
- **SAMDeflectName, SAMNoDeflectName:** used to read/write the SAM state names. For example the state may be given the name of the detector the light will fall upon.
- **biSettleDelay:** the recommended settle delay after the SAM has changed state. *BI\_select\_wavelength* will return this value if the SAM has changed state and it is the longest delay time triggered.

### 2.7.5. TLS

TLSs can be in one of three positions. TLSs have an initial position that they go to when *Bl\_initialise* is called and also maintain a list of up to 10 records that specify what position the TLS should be in at a given wavelength. If a TLS is attached to a monochromator then a call to *Bl\_select\_wavelength* will cause the TLS to go to the position corresponding to the record with the highest wavelength that is less than or equal to the wavelength specified.

**TLSCurrentPositon:** the position that the TLS is currently in (1..3).

- **TLSWI, TLSPOS:** used to access the TLS's wavelength-position records. Both attributes are indexed by the setup number (1..10). *TLSWI* refers to the wavelength and *TLSPOS* to the corresponding position.
- **biSettleDelay:** the recommended settle delay after the TLS has changed state. *BI\_select\_wavelength* will return this value if the TLS has changed state and it is the longest delay time triggered.

### 2.7.6. MVSS

The width of an MVSS can be controlled be controlled in 2 ways. The first is in a constant Bandwidth more and the second in constant width mode. In constant Bandwidth mode the desired bandwidth is set and the dll controls the slit to give the bandwidth required for the current grating. Calls to BI\_Select\_Wavelength that result in a grating change will automatically change the slits. In Constant Width mode the Slits are set to a required width in mm. The commands are:

**MVSSConstantBandwidth:** this allows for the bandwidth to be set in nm.

MVSSConstantwidth : this allows for the width to be set in mm

**MVSSSItMode**: this set which mode the slits should operate in (mm = 0) (nm = 1)

### 2.7.7. SOB

SOBs can be in one of two defined states, energised or relaxed. SOBs have an initial state that they go to when *BI\_initialise* is called.

**SOBInitialState:** the state of the SOB after calling *BI\_initialise*. 0 is relaxed and 1 is energised. The default is 0.

**SOBState:** the current state of the SOB. 0 is relaxed and 1 is energised.

**biSettleDelay:** the recommended settle delay after the SOB has changed state.

#### 2.7.8. 228, 228A, 485 and 487 ADCs

The number of samples averaged to obtain a reading is sets the integration time. For a 228A that produces a new sample every 100ms setting 10 samples per reading means that each reading takes 1s. Increasing the integration time increases the signal-to-noise ratio. Adaptive integration makes use of this by increasing the integration time as the signal weakens.

**ADCSamplesPerReading:** how many ADC samples should be averaged per measurement. Indexed 0-6, where 0 is used for non-adaptive integration and 1-6 corresponds to current amplifier range (minimum-maximum gain, i.e. strong-weak signal).

ADCAdaptiveIntegration: 0 and 1 turn adaptive integration off and on respectively.

ADCSamplePeriod: the number of milliseconds between samples. Read only.

**ADCVolts:** this will return a reading from the ADC in volts. One of few tokens which talk directly to hardware.

**ADCAuxInput:** Used to set a 487 to the Auxiliary Input.

**ADCAuxOffset:** by default, no offset is set for the Aux Input so this must be manually taken and set.

ADCAuxVolts: this will take a reading of the 487 Auxiliary Input and return a value in volts.

#### 2.7.9. 225, 265, 267, 277, 485,487 and 477 Amplifiers

The 225, 267, 277, 485,487 and 477 amplifiers can all be pre-programmed with multiple setups. The 225 and 485 have four set-ups, corresponding to four different detectors, and the 267, 277, 487 and 477 have two set-ups.

The 267/277/487/477 set-ups consist of minimum, maximum and start ranges, input channel and the wavelength from which the set-up should be used. With the 267 and 277 the set-ups are mainly used to change channel at a given wavelength.

The 225 and 485 set-ups consists of minimum, maximum and target ranges, time constant, phase quadrant and variable, frequency mode, input channel and the wavelength from which the set-up should be used. With the 225 and 485 each set-up corresponds to a different detector, and switching between set-ups provides a quick and easy way of configuring the 225 and 485 for all of the different detectors that may be encountered in a system.

When an amplifier is in the active group a call to *Bl\_initialise* will cause the amplifier to configure itself using the current setup. In addition *Bl\_select\_wavelength* will cause the amplifier to use the set-up with the highest wavelength that is less than or equal to the wavelength specified. Amplifier settings may also be 'manually' changed via *Bl\_set*. This clears AmpCurrentSetup.

AmpGain: the gain for a specific range. Indexed by the range number (1-6 for 265, 267 and 277, 1-7 for 225), where 1 is the least sensitive range. The default values for ranges from 1 up are:

265: 10<sup>4</sup>, 10<sup>3</sup>, 10<sup>2</sup>, 10<sup>,</sup> 1<sup>,</sup> 10<sup>-1</sup> 267: 10<sup>4</sup>, 10<sup>3</sup>, 10<sup>2</sup>, 10<sup>,</sup> 1<sup>,</sup> 10<sup>-1</sup> 277: 10<sup>7</sup>, 10<sup>6</sup>, 10<sup>5</sup>, 10<sup>4</sup>, 10<sup>3</sup>, 10<sup>2</sup> 225: 10<sup>-2</sup>, 10<sup>-3</sup>, 10<sup>-4</sup>, 10<sup>-5</sup>, 10<sup>-6</sup>, 10<sup>-7</sup>, 10<sup>-8</sup>

These values are used by the DLL when calculating results for *BI\_measurement* and should not normally have to be altered.

- **AmpChannel:** input channel for a specified set-up, indexed by the set-up number. Legal values are 1 and 2; the default is 1.
- **AmpMinimumRange:** minimum range to be used for a specified set-up, indexed by the setup number. Legal values are 1..number of ranges, where 1 is the least sensitive range. The default is 1. Setting this to the same value as AmpMaximumRange prevents autoranging.
- AmpMaximumRange: maximum range to be used for a specified set-up, indexed by the set-up number. Legal values are 1..number of ranges, where 1 is the least sensitive range. The default is 1. Setting this to the same value as AmpMinimumRange prevents autoranging.
- AmpStartRange: (265, 267 and 277 only) the start range to be used for a specified set-up, indexed by the set-up number. *Bl\_initialise* sets the amplifier to the start range for the current set-up. Legal values are 1..number of ranges, where 1 is the least sensitive range. The default is 1.
- AmpCurrentSetup: (267, 277 and 225 only) the current set-up. Legal values are 1 and 2 for 267 and 277, 1..4 for 225. This attribute is set to 0 if AmpCurrentRange or AmpCurrentChannel are accessed using Bl\_set.
- AmpCurrentRange: the current range. Legal values are 1..6 for the 265/267/277, 1..7 for the 225.

AmpCurrentChannel specifies the current input channel. Legal values are 1 and 2.

- **AmpUseSetup:** (267, 277, 225 only) the wavelength from which the set-up for detector n should be used, indexed by set-up number. If more than one set-up is specified for a particular wavelength then the one with the lowest value of n is used. The default is 0 (0nm).
- **biSettleDelay:** the recommended settle delay after the amplifier has changed one of its settings. *Bl\_select\_wavelength* will return this value if the amplifier has changed range or set-up and it is the longest delay time triggered.

The following apply only to the 225:

- A225TargetRange: the target range to be used for a specified set-up, indexed by the set-up number. Legal values are 1..7, where 1 is the least sensitive range. The default is 1.
- **A225PhaseQuadrant:** the phase quadrant to be used for a specified set-up, indexed by the set-up number. Legal values are 1..4, corresponding to phase quadrants of 0°, 90°, 180° and 270°. The default is 1 (0°).
- **A225PhaseVariable:** the phase variable to be used for a specified set-up, indexed by the setup number. The phase variable is a real number in the range 0..102.4. The default is 0.
- A225TimeConstant: the time constant to be used for a specified set-up, indexed by the setup number. Legal values are 1..7, corresponding to time constants of 10ms, 30ms, 0.1s, 0.3s, 1s, 3s and 10s. The default is 1 (10ms).
- A225fMode: the frequency mode to be used for a specified set-up, indexed by the set-up number. Legal values are 1 and 2, corresponding to frequency modes f and 2f. The default is 1 (f).

### 2.7.10. Motorised Stages

A Motorised Stage is simply a single motor, often used for positioning in goniometer systems.

- **MotorPosition: used to either get the current position, or tell the motor to move to a new** position. Value given is in steps.
- **MotorStop:** when MotorPosition attribute is set and passed an index of 1, the motor is not polled to work out when the motor has finished moving, giving the user an opportunity to call MotorStop if the motor needs to stop prematurely.

### 2.7.11. EBox Monitor

**EboxReadHv:** The High Voltage value, converted to volts using the gradient HV, target HV and counts at Target HV attributes.

**EboxReadTemp:** The Temperature value, converted to degrees celsius using the gradient temp, target temp and counts at target temp attributes.

**EboxReadHvRaw:** A single unconverted HV reading in ADC counts.

EboxReadTempRaw: A single unconverted temperature reading in ADC counts.

**EboxWait:** The waiting period between repeated readings of the HV or temperature values.

EboxRepeats: The number of repeats to perform for HV and temperature readings.

EboxCountsAtTargetTemp: The ADC counts at the given target temperature.

EboxGradientTemp: The gradient of the ADC counts vs temperature line.

**EboxTargetTemp:** The target temperature value in degrees celsius.

**EboxCountsAtTargetHv:** The ADC counts at the given target high voltage.

**EboxGradientHv:** The gradient of the ADC counts vs high voltage line.

**EboxTargetHv:** The target high voltage value in volts.

### 2.7.12. Camera Attributes

**CameraAutoRange:** Whether the Camera automatically changes integration time based on the maximum value across the array during a measurement.

- **CameraIntegrationTime:** The Camera integration time in milliseconds. Can be retrieved or set as appropriate.
- **CameraBeta:** The Camera Beta. This value is used along with Phi to calculate the wavelength scale across the camera array.
- **CameraPhi:** The Camera Phi. This value is used along with Beta to calculate the wavelength scale across the camera array.

### 2.7.13. System Attributes

These attributes do not belong to any particular components, but control how the system model takes measurements. In particular these attributes can be used to set-up adaptive integration in DC systems. This can speed up measurements by allowing a short integration time to be used for high signals, where the signal-to-noise ratio is already very large, and longer integration times for low signals, where the signal-to-noise ratio needs to be improved.

- **SysStopCount:** the autozero stop-count value for AC systems. When the zero offset is being calculated for an AC system (*BI\_zero\_calibration*) the DLL samples the ADC until the difference between two consecutive readings is equal to or less than this value. This ensures that the system has settled. The default is 1.
- **SysDarkIIntegrationTime:** the integration time (s) for the dark current reading in DC systems (AC systems have no dark current). The default is 5 (5s); this is equivalent to 50 samples with a 228A. Calculating the dark current is part of the zero calibration routine (*BI\_zero\_calibration*).
- **Sys225\_277Input:** the input on the 225 that the output from the 277 is connected to. The DLL uses this information to work out which combination of amplifiers is being used when it calculates the result for *BI\_measurement*. The default is 1 (i.e. 277 output to 225 input 1).

### 2.8. USING THE DLL

This section gives hints and tips on using the DLL with different programming languages and packages. As a rough guide, an application would normally call the DLL functions in the following order:

#### DLL start-up and hardware initialisation

- i) BI\_build\_system\_model
- ii) Bl\_load\_setup
- iii) BI\_initialise
- iv) BI\_park
- v) BI\_get / BI\_set (if required)
- Measurement cycle
- vi) BI\_zero\_calibration
- vii) BI\_select\_wavelength (pause using returned settle delay)
- viii) BI\_automeasure
- ix) return to step vii)

When a DLL is loaded into memory it does not have its own stack for function calls and local variables but uses the client application's stack instead. Although the spectroradiometer control DLL is not particularly stack-hungry this should be noted in case any stack errors do occur. In normal use a client stack size of 8Kb has been found sufficient (depending on how much the client uses, of course).

The serial port used by the DLL to communicate with the hardware is set in the system configuration file (hardware\system.cfg on the SDK disk). If you wish to use a different serial port then you will need to edit the line

TAS016 comms port Comn

so that COMn refers to the port that you wish to use.

The file bendll.h contains declarations for the functions in the DLL. Any source files that use functions from the DLL should include the following line:

```
#include "bendll.h"
```

The file dlltoken.h is the C/C++ token definition file. Any source files that use  $BI_get$  or  $BI_set$  need to know about attribute tokens. This is done by making sure that they contain the line:

```
#include "dlltoken.h"
```

The file dllerror.h is the C/C++ error code definition file. Any source files that need to know about error codes should contain the line:

```
#include "dllerror.h"
```

When creating a module definition file for your application a stack size of no less than 8Kb should be specified otherwise you may encounter stack errors.

The same as for C (see previous section) but the macro \_\_\_\_\_cplusplus must be defined to ensure that the DLL functions are declared correctly (some compilers do this automatically).

The file dlltoken.pas is the Pascal attribute token definition file. This should be included in any source files that use BI\_set or BI\_get (please refer to your compiler manual for instructions on how to do this). The file dllerror.pas is the Pascal error code definition file. This should be included in any files that need to know about error codes. Your compiler should also be told to allocate at least 8Kb to the stack or you may encounter stack errors.

The Visual Basic attribute token definition file and error code definition file are dlltoken.bas and dllerror.bas respectively; these need to be included in any project using the DLL. In addition the file bendll.bas contains declarations of all of the DLL functions. Including this as a

code module in a project ensures that the DLL functions are available to all other code and form modules in the project.

#### 3.1. WINDOWS ERRORS

These are error messages that Windows may report when it encounters a problem in running an application using the spectroradiometer control DLL. They appear in dialogue boxes that lock the system until dismissed.

#### File Error. Unable to find <name>.DLL.

Ensure that <*name*>.dll is correctly installed.

#### Error. Runtime error 202 at xxxx:xxxx.

This is a stack overflow error; for details on DLLs and the stack see Using The DLL.

#### Error. Runtime error 006 at xxxx:xxxx.

This error can occur when an application using the DLL has crashed and is reexecuted. When an application crashes Windows does not unload the DLL and it is left in memory in an undefined state. Trying to then access the DLL again will almost always fail, producing this error. If your application crashes the safest course of action is to restart Windows; unfortunately there is no other way of removing the DLL from memory.

### 3.2. HARDWARE CONTROL PROBLEMS

#### All of the DLL functions return BI\_error and the hardware will not respond.

Ensure that you have called BI\_build\_system\_model and that it has succeeded (i.e. returned BI\_OK).

#### BI\_initialise keeps failing

Ensure that you have called BI\_build\_system\_model and that it has succeeded (i.e. returned BI\_OK). For IEEE TM300/DTM300 systems ensure that all cable connections are secure and that the hardware is switched on.

#### The monochromator does not go to the correct wavelength.

For TM300/DTM300 systems check:

- •You have called BI\_initialise and it has succeeded,
- •You have called BI\_park and it has succeeded,
- •The DLL has the correct z-ord and  $\alpha$  values for your monochromator.

For non-parking M300/DM150 systems check: •The DLL has the correct dial reading.

For self parking M300/DM150 systems check: •You have called BI\_park and it has succeeded.

#### The filter wheel does not go to the correct position.

For TM300/DTM300 systems check:

- •You have called BI\_initialise and it has succeeded,
- •You have called BI\_park and it has succeeded.

For non-parking M300/DM150 systems check: •The DLL has the correct filter position.

For self-parking M300/DM150 systems check: •You have called BI\_park and it has succeeded, For all systems check:

•The DLL has the correct filter value for each position (any empty positions should be set to a filter value of 0).

#### Bl\_autorange keeps failing.

BI\_measurement will fail if there is a problem reading the ADC. This is often caused by a communications error or a sudden change in the input to the ADC.

•An intermittent communications error can be caused by a bad connection on the IEEE bus; check the cable connections.

•The ADC can experience a sudden change in input when the amplifier(s) is(are) ranged or the filter wheel changes position. Ensure that all delays are long enough and are actually used.

#### Bl\_zero\_calibration keeps failing.

Bl\_zero\_calibration will fail for the same reasons as Bl\_autorange; refer to the previous answer.

#### BI\_measurement keeps failing.

BI\_measurement will fail for the same reasons as BI\_autorange; refer to the previous answer.