DH-Si Silicon detector (incl. -cal, -QC and -LP options) User Notes



Overview

The Bentham DH-Si houses a 10x10mm active area silicon detector, which may either mounted to the exit slit of Bentham monochromators, or used free- standing.

This device is responsive over the spectral range 200-1100nm.

The DH-Si can be operated in either the DC or AC regimes; in the latter instance, a chopping frequency of 225Hz is suggested.

Should a calibration be purchased (_cal), please see the accompanying certificate for calibration results.

The silicon photodiode here employed is UV enhanced, and is fitted, as standard, with a quartz window. An optional resin window is available for enhanced NIR (yet compromised UV) response.

Electrical

This device is operated in the photovoltaic mode, without bias, and may be input into any Bentham current amplifier/ ORM via BNC. The pin configuration is inner pin cathode, outer sheath anode.


Mechanical

The standard DH-Si housing has hole patterns for attachment to the slits of any Bentham monochromator. The QC variant is adapted to the slit via a dove-tailed plated for easy removal of the detector. The LP housing is generally used in free space. When mounting the device to a monochromator exit slit, the fitting of a rubber O-ring around the detector should be ensured to prevent light leakage.



The device specifications are as follow:-

Typical Responsivity of Standard dh_Si



Typical detector responsivity

Active area:	10 x 10 mm
Peak sensitivity wavelength (typ.):	960 nm
Peak Responsivity (typ.)	0.52 A W ⁻¹
Dark current I _D (typ.)	2 pA
Photosensitivity vs. temperature	up to 1%/ °C at band edge
Max operating temperature:	-40 to +70 °C *
Max shunt resistance (typ.):	0.2 GΩ
NEP:	1.8 x 10 ⁻¹⁴ W/Hz ^½
Window material:	Quartz (Resin Optional)

WEEE statement:

Bentham are fully WEEE compliant, registration number is WEE/CB0003ZR. Should you need to dispose of our equipment please telephone 0113 385 4352 or 4356, quoting account number 135419.



Bentham Instruments Limited 2 Boulton Road, Reading, Berkshire, RG2 0NH, UK Tel: +44 (0)118 975 1355 Fax: +44 (0)118 931 2971 Email: sales@bentham.co.uk Internet: www.bentham.co.uk